



Bio Molecular Systems Pty. Ltd.
ABN 55141467754
5 – 7 Tonka Street
Yatala QLD 4207
AUSTRALIA

August 12, 2024

Re: Calibration and Warranty Certificate

To Whom It May Concern,

The Mic has an aluminium rotor, round and conical, spinning in an induction field to achieve heating; and cooling is achieved with a fan forced air stream. The rotor is surrounded by the induction coil, so when activated only the rotor gets hot. Unlike block cyclers that drift over time the Mic has NO DRIFT. Therefore, it should not require temperature calibration, as long as the instrument is used in a manner according to the recommendations set out in the user manual.

With a fixed optical path and no moving parts there is never any optical alignment, crosstalk compensation or optical calibration required. It is for these reason that CALIBRATION IS NOT REQUIRED and we are able to provide 3 years factory warranty. This warranty covers any defects in material and workmanship with the exceptions stated below.

How Long Coverage Lasts

This warranty runs for THREE (3) YEARS from the warranty start date ("Warranty Period").

The warranty start date is from the date of the customer invoice for any instrument that proves to be defective within the coverage of the Warranty Period. At BMS's discretion, the replacement instrument will be either a new or refurbished instrument. The defective instrument must be returned within 14 days of receiving the replacement instrument, or the replacement instrument will be invoiced at standard list price. Bio Molecular Systems will cover all international freight costs.

Extended Warranty

The standard warranty may be extended by an additional 1 or 2 years, to give a 4 or 5-year warranty respectively. Extended warranty must be taken out within 90 days of the warranty start date.

To Make a Claim

All claims must be arranged through the distributor for your region.

Kind Regards,

Valin Reja

Dr. Valin Reja
Chief Scientist | Head of Sales and Marketing